

## PhotonMAX: 512B



The PhotonMAX: 512B from Princeton Instruments/Acton is the most advanced, no-compromise electron-multiplying CCD (EMCCD) camera on the market today. Not only is it capable of single-photon sensitivity at high frame rates, it minimizes dark current by thermoelectrically cooling a back-illuminated detector up to -80°C without any water assistance. Furthermore, its all-metal, hermetic vacuum seals are guaranteed for the lifetime of the camera- only such guarantee in the industry. Its ultra-precise electronics are designed to deliver the most stable results for the most demanding quantitative applications. A 512 x 512-pixel, frame-transfer EMCCD and 10-MHz, 16-bit digitizer provide the performance needed to measure fast kinetics. In fact, the PhotonMAX is a true "2-in-1" imaging solution, boasting dual amplifiers that make the camera ideal for applications requiring either high speed or long exposures. The state-of-the-art camera interface offers circular buffers for real-time frame access and focus. It also offers optional fiber optic data interface for remote operation.

**Applications:** Single-molecule detection, spectroscopy, chemi-luminescence, astronomy, adaptive optics, hyperspectral imaging, phosphor imaging and tomography

Features	Benefits
<b>On-chip multiplication gain</b>	Low-noise, impact-ionization process for single-photon sensitivity
<b>Back-illuminated CCD</b>	>90% peak quantum efficiency offers the highest available sensitivity
<b>Frame-transfer architecture</b>	No need for mechanical shutter Allows 100% duty cycle imaging for particle-tracking applications
<b>Deep cooling</b>	Thermoelectric air cooling to -80°C minimizes dark current and allows long exposure times No need for bulky chilled-water circulators that cause condensation around the experiment optics For vibration sensitive environments, a compact room temperature liquid circulator is available as an option
<b>Lifetime vacuum</b>	Permanent, all-metal vacuum seals guaranteed for lifetime of the camera Maintenance-free operation
<b>Single optical window</b>	Vacuum window is the only optical surface between incident light and the CCD surface No losses due to multiple optical surfaces
<b>Dual amplifiers</b>	Individually optimized signal chains for a true 2-in-1 camera configuration Use the camera for applications that require high speed ("on-chip multiplication gain" amplifier) or long integration ("traditional" amplifier)
<b>16-bit digitization</b>	Wide dynamic range to capture dim and bright signals in a single image Offered at all speeds. Adjustable Analog gains to access full well capacity of the CCD.
<b>10- and 5-MHz readout</b>	Video rates at full-frame resolution (use ROI for hundreds of frames per second)
<b>1-MHz readout</b>	Slower speed of the "traditional" amplifier yields the performance of a traditional slow-scan CCD camera
<b>Kinetics readout mode</b>	Powerful readout mode offers microsecond time resolution between sub-frames
<b>PCI interface</b>	The best data-interface design in the industry
<b>Fiber optic interface (optional)</b>	Optional fiberoptic data interface for distances up to 300 meters, ideal for hazardous environments Compatible with Windows® 2000/XP, Linux®, and Mac OS X
<b>Software interface</b>	Universal interface for easy custom programming Real-time focus and image access via circular buffers Full-feature LabVIEW™ VIs available
<b>Triggering</b>	Built-in, multiple-trigger modes No external trigger boxes necessary Acquire images using single trigger, multiple triggers, or bulb trigger (allows control of exposure time via external TTL level)
<b>I/O signals</b>	Expose, readout, and shutter signals for synchronization with external devices Programmable input/output ports for complete control
<b>C-mount</b>	Easily attaches to microscopes, standard lenses, or optical equipment


## PhotonMAX: 512B Specifications


<b>CCD image sensor</b>	e2v CCD97; back-illuminated, frame-transfer CCD with on-chip multiplication gain	
<b>CCD format</b>	512 x 512 imaging pixels 16 x 16- $\mu$ m pixels 8.2 x 8.2-mm imaging area (optically centered)	
	<b>"On-chip multiplication gain" amplifier</b>	<b>"Traditional" amplifier</b>
<b>Read noise (typical)</b>	~40 e- rms @ 5 MHz ~55 e- rms @ 10 MHz Read noise effectively reduced to <1 e- rms with on-chip multiplication gain enabled	8 e- rms @ 1 MHz 15 e- rms @ 5 MHz
<b>Full well (typical)</b>	800 ke-	200 ke-
<b>Non-Linearity</b>	<2%	<1%
<b>Analog gain (typical)</b>	12, 6, 3 e-/p/s	4, 2, 1 e-/p/s
<b>Deepest cooling temperature</b>	-70°C $\pm$ 0.05°C	-80°C $\pm$ 0.05°C
<b>Dark current@ -70°C</b>	0.005 e-/p/sec (typical) 0.01 e-/p/sec (maximum)	
<b>Spurious events (typical)</b>	0.005 e-/pixel/frame clock induced events measured with 33msec exposure time and ~1000x multiplication gain	
<b>On-chip multiplication gain</b>	1 to 1000x Software controlled in 4,096 steps	
<b>Digitization</b>	16 bits @ 10 MHz, 5 MHz, and 1 MHz	
<b>Vertical shift rate</b>	2 $\mu$ sec/row*	
<b>Binning</b>	Flexible binning capabilities in parallel direction; 1 through 6 binning in serial direction	
<b>Operating environment</b>	0 to 30°C ambient, 0 to 80% relative humidity noncondensing	

Notes: All specifications subject to change.

\*All performance parameters such as full well, charge transfer efficiency, read noise and image quality are optimized at 2 usec/row. Contact factory for faster vertical speed operation. It must be noted that faster vertical shift speeds adversely affect the quantitative performance.

## QE Curve





Frame Rate

Binning	Region			
	512 x 512	256 x 256	128 x 128	64 x 64
1 x 1	29	54	95	155
2 x 2	56	95	155	227
4 x 4	98	155	227	295
6 x 6	130	195	262	329

(Frames per second)

Note: Frame rates are measured at 10 MHz digitization and 2usec/row vertical shift speed

Princeton Instruments



www.piaction.com

email: moreinfo@piaction.com

USA +1.877.4 PIACTION | Benelux +31 (347) 324989

France +33 (1) 60.86.03.65 | Germany +49 (0) 89.660.779.3

Japan +81.3.5639.2741 | UK +44 (0) 28.38310171